

# Beltone Commence™

Introducing  
hearing clarity and  
connectivity for all



For professionals

# Beltone Commence™

Introducing  
hearing clarity and  
connectivity for all



For professionals

# Clarity and ease

Offer even more clients the proven all-around sound quality<sup>1-4</sup> of Beltone Commence – available as an extended lineup with RIEs, including the miniRIE, and now Customs.



## Best-in-class solutions

Our small, all-weatherproof\* designs give all-day wearing comfort and all-day battery power<sup>5</sup> for a great sound experience with leading feedback management.<sup>6</sup> You can offer rechargeable and battery-type styles.

### Smart Gain

For clear and comfortable hearing that automatically adapts to rapidly changing sound environments<sup>2</sup>

### CrossLink Directionality

For great hearing clarity all day<sup>1</sup>

### Impulse Noise Reduction



For greater listening comfort<sup>4</sup>

### Effortless connections

Your clients can enjoy best-in-class Bluetooth® connectivity with Bluetooth Low Energy (LE) Audio and connect to Auracast™ broadcast audio.\*\* Even the new ITC Custom hearing aid offers Auracast!

### Auracast Assistant

Find and select Auracast broadcasts with a tap using the Auracast Assistant in the Beltone HearMax™ app.



New Customs

### New family of Customs

Experience our biggest upgrade of sound quality, wearing comfort and connectivity since our previous Essential generation.

### Custom-made options at your fingertips

Recommend the CIC hearing aid for an easy wear & forget fit, or the wireless ITC design for a comfortable all-day fit and a connected experience. The ITC hearing aid also features an updated generation of the microchip to enable a richer sound.



### Introducing AllFit domes


Explore our range of AllFit domes in 19 styles and sizes, designed for the best comfort and acoustic performance, and to fit at least 97% of ear types and sizes.<sup>7</sup>

### Elevated fitting experiences

Our flex demo devices are made for showcasing the sound features – and for loaning out, so your clients are never without hearing aids.\*\*\* Also, exchangeable mic filters on all our styles\*\* give greater flexibility and reduce time spent sending devices in for servicing.

1. Groth (2016), 2. Groth & Cui (2023), 3.Haastrup (2012), 4. Sjolander et al. (2019), 5. All-day battery covers a minimum of 12 hours of use based on average use data for RIE styles. Battery drain measured with streaming and in hearing-in-noise program. GN Proprietary data on file (2025), 6. Compared to leading solutions. Cui & Groth (2025) 7. GN Proprietary data on file (2025), \*Fully IP68-rated, \*\*Not for CIC (COM4,3,2), \*\*\*Only for RIEs

# Styles and colours



## Key feature overview

For a detailed feature overview across styles, please refer to the data sheets

Features	COM4	COM3	COM2
<b>Sound Quality</b>			
WDRC bands	★★★★☆ 12	★★★★☆ 12	★★★★☆ 10
Environment recognition	★★★★☆ (Smart Gain)		
Directionality	★★★★☆ (CrossLink Directionality 2)*	★★★★☆ (Speech Spotter)*	★★★★☆ (Speech Spotter)*
Speech intelligibility in noise	★★★★☆ (Smart Beam Steering)*	★★★★☆ (Smart Beam Steering)*	★★★★☆ (Adaptive Directionality)**
Noise reduction	★★★★☆ (2 settings)	★★★★☆ (1 setting)	★★★★☆ (1 setting)
Comfort	★★★★☆ (2 settings)	★★★★☆ (1 setting)	★★★★☆ (1 setting)
Feedback management	★★★★★ (Feedback Eraser 2 w/Music Mode)	★★★★★ (Feedback Eraser 2)	★★★★★ (Feedback Eraser 2)
Low frequency boost	★★★★★ (3 settings)**	★★★★★ (3 settings)**	★★★★★ (3 settings)**
Sound shifter	Yes		
<b>Connectivity</b>			
Hands-free calls	★★★★★ Direct streaming for iPhone and iPad** and Hands-free calls*. Connecting hands-free to compatible Bluetooth 5.3 or later devices, including select Android phones.		
Public connectivity	★★★★★ Supports Auracast for new broadcast opportunities.** Telecoil built in RIE64 (and in Multi-Mic+).		
Accessories	★★★★★ Multi-Mic+ and TV Streamer+, Remote Control, TV Streamer 2, Beltone HearMax App.**		
<b>Fitting Flexibility</b>			
Flexible demo devices	Yes		
Gain handles	12	12	10
Programs	4	4	4
Patient experience profiles	Yes		

**Worldwide Headquarters**  
GN Hearing A/S  
Lautrupbjerg 7  
DK-2750 Ballerup  
Denmark  
Tel.: + 45 45 75 11 11  
gnhearing.com

**Australia**  
GN Hearing Australia Pty Ltd  
Suite 1001, Level 1  
1 Julius Avenue  
North Ryde, NSW 2113  
Australia  
Tel.: (free) 1800 658 955  
gnhearing.com

\*Not for ITC & CIC  
\*\*Not for CIC

Learn more at  
gnhearing.com



© 2026 GN Hearing A/S. All rights reserved. Beltone is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone and iPad are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google LLC. The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast word mark and logos are trademarks owned by Bluetooth SIG, Inc.